

ROBOLEX



POWER	AC110 ~ 220V, 50/60Hz
POWER CONSUMPTION	MAX 500VA
RF OUPUT(MAX)	20W(at 500Ω)
RF OUTPUT FREQUENCY	1MHz ~ 1.8MHz(RV95, 60) 2MHz ~ 3MHz(RV25) (Automatic frequency adjustment)
RF OUTPUT TYPE	MULTI(RV95, 60), BI(RV25)
VACUUM	MAX 0.8bar(600mmHg) (AC Pump Cont. Max 2.1Bar)
LASER	655nm(3mW)-RV95,60
ULTRASOUND FREQUENCY (Optional)	38Khz

ROBOLEX







ROBOLEX

ROBOLEX is an effective and space saving body contouring device combining Vacuum suction, Ultrasound Cavitation, Low Level Laser and Multi/Bi-Polar Radio Frequencies.

Lipocyte destruction by using Ultrasound Cavitation.

Dissolving destroyed lipocyte by using Multi/Bi-Polar Radio Frequency and cellulite reduction by Low Level Laser.

Drain the dissolved lipocyte rapidly from the body to outside by the lymph drainage and blood circulation and maintain the elasticity on treated area with Vacuum.



Dissolve fat where under 5mm from the skin by 655nm laser RF infiltrates into $5\sim15$ mm from the skin and spread oxygen into intracellular by heating Soften skin with negative pressure made by vacuum and help to distribute energy evenly

APPLICATION _____



Skin's electricity recover, Lifting

Skin Tightening

Wrinkle Improvement

Pore Improvement

ROBOLEX



Body contouring

Management of Patient who took liposuction

SUMMARY -

METHOD

Obesity care by Multi-functional treatment (RF + Endomologie + Laser)
Optional item : Ultlasound cavitation.

FEATURE

Completion of wanted body shape by selective reduction of fat cell and Cellulite with non invasive treatment.

THEORY

Destruction of lipocyte by combination operation with negative pressure in lipocyte & ultrasound cavitation.

SAFETY

Complementary lipolysis treatment way in safe to prevent the risk by liposuction and to get rid of uncertainty by injection therapy.

EFFECT

Confirmation of reduction for size & body shape by 3~4 times treatment and it preveats the wrinked line on the reduced part after treatment.

